CM90

Directly Compressible granulated Calcium Carbonate with Maltodextrin for supplement tablets



With our special granulation technology, we provide our mineral supplement customers with a highly compressible high assay (91%) calcium carbonate.

SPI Pharma

An ABF Ingredients Company

Ideally CM90 is used for swallow tablets due to its high density and compressibility. However, it can also be used for chewable, gummies and other tablet forms.



Why use CM90 for swallow supplement tablet?

- High density DC Calcium Carbonate is highly compressible, allowing very small easier to swallow tablets
- Excellent flow
- CM90 has excellent compactability and binding properties allowing the formulator to reduce the amount of other binders in the formula; making room for additional actives or simply making a smaller tablet

CM90

Excellent for tableting by direct compression, saving time and money

Excellent compactibility and compressibility

• At compression force between 10-30 kN, the CM90 tablets have ideal hardness and thickness for swallow tablets

Good flowability

- Our granulation technology provides consistent PSD from batch to batch and that ensures good flowability and good tableting
- Consistent feed to the tableting machine





Versatility

- CM90 is a versatile ingredient that can be used as an active or as an excipient. It provides buffering properties (in a pH 9 environment)
- Can be used for direct compression and wet granulation
- Can be used in many different final application forms - swallow, chewable, sachet, gummies, toffees, and effervescent forms with smooth fizzing



ingredients in the formulation.

Starting Formula using CM90 in a swallow Calcium and Vitamin tablet

	LABEL CLAIN	1 per Ta
Calcium Carbonate	2500 mg CaCo	
	equivalent to: 1000 mg of Co	
Vitamins D3	800 IU	
INGREDIENTS	% by weight	Qty
CM90 at 91% CaCO ₃	96.73 %	27
Vitamin D3 (100 IU/mg)	0.28 %	
Crospovidone	2.00 %	5
Magnesium Stearate	0.99 %	
	Total Table	et Weigl

Manufacturing Guidance

- 1. Screening / Weighing: screen through # 3 and weigh the ingredients listed in the above recipe
- 2. Blending: blend the ingredients for approx 10 minutes on a Turbula Blender Type
- 3. Tableting: compress the blend to the targe characteristics on a rotary tablet press

Elemental Impurities

Operating in a highly regulated industry, we und the importance of minimal elemental impuritie customers. We perform a thorough risk assessm our products to ensure regulatory compliance. has invested approximately \$1M of testing tech resources and external testing to ensure our cus peace of mind that our products are below the and only require annual confirmation testing. C depend on our certification and can alleviate their burden of testing.



SPI Pharma's expertise will help you increase your production rate especially for smaller tablets, by providing you with a high flow product that blends consistently and homogeneously with other

ablet
CO ₃
Calcium (Ca²+)
y per Tablet
747.25 mg
8.00 mg
56.80 mg
27.95 mg
ght 2840.00 mg
30 mesh sieve
ximately
5
eted tablet
derstand
es for our
nent on all of SPI Pharma
nology,
stomers the
threshold
ustomers can

Regulatory Compliance[†]

CM90 Calcium Carbonate		
\checkmark	Technical File	
\checkmark	Non-GMO *	
\checkmark	Gluten Free	
\checkmark	Halal (for ingredients)	
\checkmark	Kosher	
\checkmark	Aflotoxin Free	
\checkmark	cGMP	
\checkmark	Nano Material Free	
\checkmark	Components Origin	
\checkmark	Latex Free	
\checkmark	Absence of Residual Solvents	
\checkmark	Absence of Phthalates	
\checkmark	Food Compliance	
\checkmark	Heavy Metals Compliance	
\checkmark	USP	
\checkmark	BSE/TSE	
\checkmark	Absence of Allergens	

* Available for Calcium Carbonate No Non-GMO for Maltodextrin

† Documents available upon request





All information and statements given in this brochure are believed to be accurate at the time of publication. However, neither SPI Pharma nor any of their affiliates make any representations or warranty with respect thereto, including, but not limit to, any results obtained in the processing of the products by customers or any third party. All information and statements are intended for persons having the required skill and know-how and do not relieve the customer or user from verifying the suitability of information and statements given for a specific purpose prior to use of products. It is entirely the obligation of the costumer or user to comply with applicable laws and regulations, and also with all patent or other intellectual property rights of third parties.

SPI PHARMA EXPRESSLY DISCLAIMS ANY REPRESENTATIONS OR WARRANTIES OF ANY KIND. WHETHER EXPRESS OR IMPLIED, AS TO THE ACCURACY, CURRENCY, COMPLETENESS AND/OR THE MERCHANTABILTY OR FITNESS FOR A PARTICULAR PURPOSE OF ANY INFORMATION CONTAINED IN THIS BROCHURE AND/OR ANY PRODUCT DESCRIBED OR PROMOTED IN THIS BROCHURE, INCLUDING WARRANTIES WITH RESPECT TO INFRINGEMENT OF ANY PATENT, COPYRIGHT, OR OTHER RIGHTS OF A THIRD PARTY. We reserve the right to change product specification and not specified properties of the products without prior notice.

> October 2019 All rights reserved © 2019 SPI Pharma